

Fig. S1. Cytotoxicity of PYC in H292 cells. The graph bars indicate the cytotoxicity of PYC with H292 cells. Cells were seeded into 96 well plates, thus PYC (2.5-5-10-20-40 $\mu\text{g/ml}$) or AMP (2-Amino-4-methylpyridine; 5 μM) was treated. After treating the MTT solution of 5 $\mu\text{g/ml}$, plates were analyzed using a multifunctional plate reader (SPARK 10 M, Tecan, Switzerland) at OD 570 nm.

